






## Connex Vital Sign Monitor 6000 Series (CVSM) Competency Test

Name: \_\_\_\_\_ Unit: \_\_\_\_\_ Date: \_\_\_\_\_

- 1) The proper way to power down the Connex Vital Signs Monitor 6000 Series, would be to:
  - a) Do not touch the screen for 2 minutes and it will automatically power down
  - b) Touch the Battery icon located on the upper right hand section of the touch screen and select "Power Down"
  - c) Touch Review tab, then Send
  - d) None of the above.
  
- 2) True or False: The Connex Vital Signs Monitor 6000 Series is intended to be used on adult, pediatric, and neonatal patients.
  
- 3) The monitor will change into Direct temperature mode in all of these ways *except*:
  - a) 60 seconds after the thermometer probe is removed from the probe well and exposed to air.
  - b) After you complete a Predictive mode measurement, then touch  to switch from Predictive to Direct mode.
  - c) Not using a Welch Allyn thermometer probe cover.
  - d) If you have a patient whose body temperature is below the normal temperature range.
  
- 4) True or False: You can configure the monitor to take automatic NIBP measurements at variable intervals.
  
- 5) The CVSM 6000 LED light on the device handle will do which of the following during a high-priority alarm?
  - a) The handle will remain a constant amber color.
  - b) The handle will flash an amber light.
  - c) The handle will flash red.
  - d) The handle will show an alternating green and red color.
  
- 6) True or False: During an SpO<sub>2</sub> measurement, the displayed pulse rate is derived from the SpO<sub>2</sub> sensor. If SpO<sub>2</sub> is not available, the pulse rate is derived from NIBP.

- 7) The following icon represents that all alarms audible and visual alarms are off.
- a) 
  - b) 
  - c) 
  - d) 
- 8) You can manually enter patient parameters except which of the following?
- a) Height
  - b) Weight
  - c) Respirations
  - d) Age
- 9) True or False: The monitor enables the scanning of patients' and/or clinicians' barcodes to enter identification information.
- 10) True or False: Cavi-Wipes™, Sani-Cloth® Plus, 10% chlorine bleach solution, and 70% alcohol are all approved cleaning solutions for the device.



## CVSM 6000 Competency Test Answer Sheet

Name: \_\_\_\_\_

Unit: \_\_\_\_\_

Date: \_\_\_\_\_

1. \_\_\_\_\_

2. \_\_\_\_\_

3. \_\_\_\_\_

4. \_\_\_\_\_

5. \_\_\_\_\_

6. \_\_\_\_\_

7. \_\_\_\_\_

8. \_\_\_\_\_

9. \_\_\_\_\_

10. \_\_\_\_\_

## CVSM 6000 Competency Test Answer Key

1. B
2. True
3. C
4. True
5. C
6. True
7. A
8. D
9. True
10. True