

FlexiPort® EcoCuff™


A cuff designed to control cross-contamination, costs and the volume of hospital waste in a single-patient-use setting.

EcoCuff

- Reduce the Risk of Cross-Contamination
- Control Cost
- Capable of being recycled*

Learn more
on page 94



 The EcoCuff is capable of being recycled. Although the practice is not yet widespread, a number of hospitals have implemented the practice of recycling medical-use polypropylene plastics.



Sphygmomanometers

767 Wall & Mobile Aneroids.....	72
Platinum Series DS58™ Hand Aneroids	74
Platinum Series DS48 Integrated Aneroids	75
Platinum Series DS48A Pocket Aneroids	76
Gold Series DS66 Trigger Aneroids.....	77
Silver Series DS45 Integrated Aneroids	78
Bronze Series DS44 Integrated Aneroids	80

Home Blood Pressure Measurement

OSZ 5 Home Blood Pressure Systems	82
---	----

Flexiport® Blood Pressure Cuffs

FlexiPort Reusable One-Piece Blood Pressure Cuffs.....	83
FlexiPort Soft Disposable One-Piece Blood Pressure Cuffs	87
FlexiPort Vinyl Disposable One-Piece Blood Pressure Cuffs	91
FlexiPort EcoCuff™ Single-Patient-Use Blood Pressure Cuffs.....	94
FlexiPort Accessories	95
Neonatal Soft Disposable One-Piece Blood Pressure Cuffs with Luer Slip Connectors (10 pack)	97

Trimline™ Blood Pressure Cuffs

Trimline Reusable Blood Pressure Cuffs.....	98
Trimline Soft Disposable Blood Pressure Cuffs	102
Trimline Vinyl Disposable Blood Pressure Cuffs.....	107
Trimline Disposable Isolation Blood Pressure Cuffs.....	109
Trimline Soft Neonatal One-Piece Blood Pressure Cuffs.....	110
Trimline Vinyl Neonatal One-Piece Blood Pressure Cuffs with Luer Slip Connectors.....	111

Two-Piece Cuffs and Bladders

One-Tube Blood Pressure Cuffs and Neoprene Bladders	113
One-Tube Neoprene Bladders.....	113
Two-Tube Blood Pressure Cuffs and Neoprene Bladders.....	114
Two-Piece Inflation Systems	114
Two-Tube Neoprene Bladders	115

Bulbs/Valves

Inflation Bulbs and Valves.....	115
---------------------------------	-----

SPHYGMOMANOMETERS

BLOOD PRESSURE MANAGEMENT



7670-01
Wall Aneroid with
Adult Cuff

767 Wall and Mobile Aneroids

Improved ergonomics and a new modern design make the 767 Aneroid an easy-to-use and attractive blood pressure aneroid. Based on Welch Allyn's near century of experience building and designing innovative blood pressure technology, it's the smart choice for accurate and reliable readings in your facility.

- Laser-engraved dial face provides unmatched accuracy
- Recessed dial face and high-contrast red pointer increase visibility and reduce parallax errors
- Lifetime calibration warranty
- New 5-leg stand provides improved stability and maneuverability
- Simple color changes give the aneroid a new look
- Mercury-free
- Latex-free

767 Wall Aneroids

7670-01 767 Wall Aneroid with Reusable Flexiport Adult Cuff

7670-01CB 767 Wall Aneroid with Traditional Two-Piece Cuff and Bladder Inflation System

7670-01C 767 Wall Aneroid with Reusable Flexiport Child Cuff

7670-02 767 Wall Aneroid with 8 ft. Coiled Tubing (cuff not included)

767 Mobile Aneroid

7670-03 767 Mobile Aneroid with Reusable FlexiPort Adult Cuff



7670-03
Mobile Aneroid

767 Wall and Mobile Aneroids: Accessories

7670-06 767 Mobile Stand and Basket Assembly Only

7670-07 Inflation System Basket, Large, White 767 Wall and Mobile Aneroids

7670-08 Crystal, 767-Series 767 Wall and Mobile Aneroids

5086-01 Premium Small Inflation Bulb (Black) 767 Wall and Mobile Aneroids

5086-03 Premium Large Inflation Bulb (Black) 767 Wall and Mobile Aneroids

5087-01 Air Release Valve 767 Wall and Mobile Aneroid

5088-01 Air Release Valve and Small Bulb 767 Wall and Mobile Aneroids

5082-167 Fitting, Metal Male Luer Slip (10 pack) 767 Wall and Mobile Aneroids

5082-169 Fitting, Metal Female Luer Slip (10 pack) 767 Wall and Mobile Aneroids

5089-39 Coil Tubing 4 ft Length 767 Wall and Mobile Aneroids

5089-13 Coil Tubing 8 ft Length 767 Wall and Mobile Aneroids

5089-13H Coil Tubing, 8 ft Length (Hospital Pack, 12/Box)

5082-23 Large Adult Inflation System with Two-Tube Bag and Cuff

5082-22 Adult Inflation System with Two-Tube Bag

5082-21 Child Inflation System with Two-Tube Bag and Cuff

5082-26 Large Adult Inflation System with Two-Tube Bag

SPHYGMOMANOMETERS

BLOOD PRESSURE MANAGEMENT



Ergonomic Inflation Bulb

DS58
DS58 Gauge Only

Platinum Series DS58™ Hand Aneroids

Maintain accurate blood pressure readings with lower lifetime cost with gear-free DuraShock™ technology. Part of our Platinum Series, this premium gauge provides unsurpassed reliability in a traditional handheld design.

- Super shock-resistant—can withstand up to 60" drop (twice the AAMI "Shock Resistant" standard) and remain in calibration
- Jewel movement for long life
- Laser-engraved dial face for accuracy
- Premium inflation system for comfort
- Lighter weight than traditional models—won't weigh down your lab coat
- Latex-free



DS58-11
Gauge with FlexiPort® Adult Cuff



DS58-MC
Gauge with FlexiPort® Family Practice Kit



DS58-PD
Gauge with Pediatric FlexiPort® Cuff Kit

DS58 Gauge Only (NEW replacement for 5098-03)

DS58-11 Gauge with Adult Durable FlexiPort Cuff and Nylon Zipper Case

DS58-11CB Gauge with Adult Cuff and Bladder and Nylon Zipper Case

DS58-MC Gauge with Family Practice Kit (Includes: Large Adult, Adult, Small Adult Child Cuffs, and Nylon Zipper Case)

DS58-MCCB Gauge with Family Practice Kit with Cuff and Bladder (Includes: Large Adult, Adult, Child-Print Cuffs, and Nylon Zipper Case)

DS58-PD Gauge with Pediatric Cuff Kit (Includes: Child-Print, Small Child, Infant and Newborn Cuffs)

DS58-PDCB Gauge with Pediatric Kit with Cuff and Bladder (Includes: Child-Print, Small Child, Infant Cuffs, and Nylon Zipper Case)

Platinum Series DS58 Hand Aneroids: Accessories

5085-10 Small Nylon Zipper Case

5085-11 Large Nylon Zipper Case

5087-10 Replacement Knob (Includes bulb and spring)

5086-06 Inflation Bulb

Platinum Series DS48 Integrated Aneroids

Maintain accurate blood pressure readings with lower lifetime cost with gear-free DuraShock™ technology, the first major innovation in aneroid technology in 100 years. Part of our Platinum Series, this integrated-style gauge and one-tube adapter plug directly into the integrated port of a FlexiPort® cuff. In use, the gauge rests on the outer part of the patient's arm.

- Super shock-resistant—can withstand up to 60" drop (twice the AAMI® "Shock Resistant" standard) and remain in calibration
- Jewel movement for long life
- Laser-engraved dial face for accuracy
- Premium inflation system for comfort
- Lifetime calibration warranty
- Lighter weight than traditional models
- Latex-free



DS48 -11
DS48 Gauge with
Adult Cuff and Case

Platinum Series DS48 Integrated Aneroids: Accessories

DS48	Gauge Only DS48 Integrated Aneroids Model Configuration
DS48-11	Gauge with Durable One-Piece Adult Cuff and Nylon Zipper Case (NEW replacement for 5090-02)
5085-10	Small Nylon Zipper Case
1SRCL	Universal Adapter with one 24" tube and a large bulb with premium valve
1SRCS	Universal Adapter with one 24" tube and a small bulb with premium valve
5082-159	FlexiPort cuff Launder Cap

DuraShock Integrated Cuffs (Gauge not included)

DS45-DUR12	FlexiPort Reusable Blood Pressure Cuff, Large Adult
DS45-DUR11	FlexiPort Reusable Blood Pressure Cuff, Adult
DS45-DUR10	FlexiPort Reusable Blood Pressure Cuff, Small Adult
DS45-DUR09	FlexiPort Reusable Blood Pressure Cuff, Child

SPHYGMOMANOMETERS

BLOOD PRESSURE MANAGEMENT



DS48-11CB
Gauge with Adult Cuff

Platinum Series DS48A Pocket Aneroids

Maintain accurate blood pressure readings with lower lifetime cost with gear-free DuraShock™ technology, the first major innovation in blood pressure technology in 100 years. Part of our Platinum Series, this lightweight pocket-style gauge has a metal clip that attaches to the cleat on a two-tube cuff. In use, the gauge rests on the inner part of the patient's arm.

- Super shock-resistant—can withstand up to 60" drop (twice the AAMI "Shock Resistant" standard) and remain in calibration
- Jewel movement for long life
- Laser-engraved dial face for accuracy
- Premium inflation system for comfort
- Lighter weight than traditional models—won't weigh down your lab coat
- Latex-free

DS48A Gauge Only (NEW replacement for 5090-03)

DS48-11CB Gauge with Durable Two-Piece Adult Cuff and Bladder and Nylon Zipper Case (NEW replacement for 5090-02)

Platinum Series DS48A Pocket Aneroids: Accessories

5085-10 Small Nylon Zipper Case



DS48A
Gauge only

Common Cuffs and Inflation Systems

5082-21 Child Inflation System with Two-Tube Bag and Cuff

5082-22 Adult Inflation System with Two-Tube Bag and Cuff

5082-23 Large Adult Inflation System with Two-Tube Bag and Cuff

5082-45 Thigh Inflation System with Two-Tube Bag and Cuff

Reuse-08-2BV One-Piece Reusable Two-tube BP Cuff, SM Child

Reuse-09-2BV One-Piece Reusable Two-tube BP Cuff, Child

Reuse-10-2BV One-Piece Reusable Two-tube BP Cuff, Sm Adult

Reuse-11-2BV One-Piece Reusable Two-tube BP Cuff, Adult

Reuse-12-2BV One-Piece Reusable Two-tube BP Cuff, Large Adult

Reuse-13-2BV One-Piece Reusable Two-tube BP Cuff, Thigh

5088-01 Air Release Valve and Small Bulb

5086-03 Large Inflation Bulb

5086-01 Small Inflation Bulb

5087-01 Air Release Valve



5085-10
Small Nylon Zipper Case

Gold Series DS66 Trigger Aneroids

Maintain accurate blood pressure readings with lower lifetime cost with gear-free DuraShock™ technology, the first major innovation in aneroid technology in nearly 100 years. Part of our Gold Series, the DS66 Trigger features an innovative trigger-release design for easy operation.

- Shock resistant to maintain accuracy with lower lifetime cost—can withstand up to 30" drop and remain in calibration
- Soft comfort grip
- Jewel movement for long life
- Laser-engraved dial face for accuracy
- Trigger deflation valve for smooth, one-finger operation
- 15-year calibration warranty
- Lighter weight than traditional models—won't weigh down your lab coat
- Latex-free

5098-30 Gauge and Multi-Cuff Kit with FlexiPort One-Piece Child-Print, Small Adult, Adult, and Large Adult Cuffs and Zipper Case

5098-31 Gauge and Bulb only

5098-27 Gauge and FlexiPort One-Piece Adult FlexiPort Cuff and Nylon Zipper Case

5098-28 Gauge and FlexiPort One-Piece Large Adult Cuff and Nylon Zipper Case

5098-29 Gauge and FlexiPort One-Piece Child Cuff and Nylon Zipper Case

5098-29CB Gauge and FlexiPort Two-Piece Child Cuff and Bladder and Nylon Zipper Case

5098-27CB Gauge and FlexiPort Two-Piece Adult Cuff and Bladder and Nylon Zipper Case

5098-28CB Gauge and FlexiPort Two-Piece Large Adult Cuff and Bladder and Nylon Zipper Case

5098-30CB Gauge and Multi-Cuff Kit with FlexiPort Two-Piece Child-Print, Adult, and Large Adult Cuffs and Bladders and Nylon Zipper Case

Gold Series DS66 Trigger Aneroids: Accessories

5086-04 Bulb Assembly

Common Cuffs

Reuse-09-1TP Child-Print [size 9]

Reuse-10-1TP Small Adult [size 10]

Reuse-11-1TP Adult [size 11]

Reuse-12-1TP Large Adult [size 12]

Reuse-13-1TP Thigh [size 13]

5082-62 Child-Print Blood Pressure Cuff and One-Tube Bladder

5082-43 Adult Blood Pressure Cuff and One-Tube Bladder

5082-44 Large Adult Blood Pressure Cuff and One-Tube Bladder



5098-31
Trigger Aneroid
DS66



5098-30
Gauge and Multi-FlexiPort Cuff Kit, Small Adult, Adult and Large Adult FlexiPort® Cuffs and Nylon Zipper Case



5098-27
Trigger Aneroid with Adult FlexiPort® Cuff and Nylon Zipper Case

SPHYGMOMANOMETERS

BLOOD PRESSURE MANAGEMENT



DS45-11
Gauge with Adult Cuff

Silver Series DS45 Integrated Aneroids

Gear-free DuraShock™ design is much more durable than traditional aneroids with gears. Part of our Silver Series, the DS45 Integrated features a unique integrated, cuff-mounted design that is compact and lightweight for ease of use.

- Shock resistant to maintain accurate blood pressure readings with lower lifetime costs—can withstand up to 30" drop and remain in calibration
- Laser-engraved dial face for accuracy
- Gauge rotates 360° for easy viewing and snaps directly into cuff for quick cuff change
- 10-year calibration warranty
- Large inflation bulb and high-quality air-release valve make inflation easy
- Latex-free

DS45 Gauge Only

DS45A Gauge with Integrated Adapter for use with Two-Tube Cuffs

DS45-09 Gauge with Durable One-Piece Child-Print

DS45-09CB Gauge with Durable Two-Piece Child Cuff and Bladder

DS45-09PCB Gauge with Durable Two-Piece Child-Print Cuff and Bladder

DS45-10 Gauge with Durable One-Piece Small Adult Cuff

DS45-11 Gauge with Durable One-Piece Adult Cuff

DS45-11C Gauge with Durable One-Piece Adult Cuff and Case

DS45-11CB Gauge with Durable Two-Piece Adult Cuff and Bladder

DS45-12 Gauge with Durable One-Piece Large Adult Cuff

DS45-12CB Gauge with Durable Two-Piece Large Adult Cuff and Bladder

Silver Series DS45 Integrated Aneroids: Accessories Common Cuffs and Inflation Systems

DS45-DUR09F	Child FlexiPort Durable One-Piece Cuff Inflation System
DS45-DUR10F	Small Adult FlexiPort Durable One-Piece Cuff Inflation System
DS45-DUR11F	Adult FlexiPort Durable One-Piece Cuff Inflation System
DS45-DUR12F	Large Adult FlexiPort Durable One-Piece Cuff Inflation System
DS45-DUR12	Large Adult FlexiPort Durable One-Piece Cuff for use with DS45 and DS48 Integrated Gauges
DS45-DUR11	Adult FlexiPort Durable One-Piece cuff for use with DS45 and DS48 Integrated Gauges
DS45-DUR10	Small Adult FlexiPort Durable One-Piece Cuff for use with DS45 and DS48 Integrated Gauges
DS45-DUR09	Child FlexiPort Durable One-Piece Cuff for use with DS45 and DS48 Integrated Gauges
5082-21	Child Inflation System with Two-Tube Bag and Cuff
5082-22	Adult Inflation System with Two-Tube Bag and Cuff
5082-23	Large Adult Inflation System with Two-Tube Bag and Cuff

Miscellaneous Accessories

5082-159	FlexiPort Cuff Launder Cap
5082-255	Calibration Connector
5085-07	Polyester Zipper Case
5086-01	Premium Small Inflation Bulb (black)
5086-03	Premium Large Inflation Bulb (black)
5087-01	Air Release Valve
5082-257	Integrated Gauge Adapter
1DSAL	DuraShock Adapter with one 24" tube and a large bulb with premium valve
1DSAS	DuraShock Adapter with one 24" tube and a small bulb with premium valve



5082-257
Integrated Gauge Adapter



5086-03
Premium Large Inflation
Bulb (black)



5087-01
Air Release Valve

SPHYGMOMANOMETERS

BLOOD PRESSURE MANAGEMENT



DS44-MC
DuraShock Multi-FlexiPort
Cuff Kit

Bronze Series DS44 Integrated Aneroids

Gear-free DuraShock™ design is much more durable than traditional aneroids with gears. Part of our Bronze Series, the DS44 Integrated features a FlexiPort® cuff-mounted design that is compact and lightweight for ease of use.

- Integrated design
- Premium inflation system
- Five-year calibration warranty
- Latex-free

DS44 Gauge Only

DS44A Gauge with Integrated Adapter for use with Two-Tube Cuffs

DS44-09 Gauge with Durable One-Piece Child Cuff

DS44-10 Gauge with Durable One-Piece Small Adult Cuff

DS44-11 Gauge with Durable One-Piece Adult Cuff

DS44-12 Gauge with Durable One-Piece Large Adult Cuff

DS44-11C Gauge with Durable One-Piece Adult Cuff and Case

DS44-MC Gauge Multi-Pack with Large Adult, Adult, Small Adult and Child Cuffs

DS44-09CB Gauge with Durable Two-Piece Child Cuff and Bladder

DS44-11CB Gauge with Durable Two-Piece Adult Cuff and Bladder

DS44-12CB Gauge with Durable Two-Piece Large Adult Cuff and Bladder

DS44-13CB Gauge with Durable Two-Piece Thigh Cuff and Bladder

DS44-11V Gauge with Durable One-Piece Adult Cuff, Includes Color Bumper Variety Pack

Bronze Series DS44 Integrated Aneroids: Accessories

DS44-DUR12	Large Adult FlexiPort, Durable One-Piece Cuff for use with DS44 Integrated Gauges
DS44-DUR11	Adult FlexiPort, Durable One-Piece Cuff for use with DS44 Integrated Gauges
DS44-DUR10	Small Adult FlexiPort, Durable One-Piece Cuff for use with DS44 Integrated Gauges
DS44-DUR09	Child FlexiPort, Durable One-Piece Cuff for use with DS44 Integrated Gauges
5082-159	FlexiPort cuff Launder Cap
5082-255	Calibration Connector
5085-09	Storage Case, Vinyl
5086-08	Bulb, Standard, Gray, Medium
5087-16	Air Release Valve, Standard
5082-257	Integrated Gauge Adapter
5082-21	Child, Inflation System with Two-Tube Bag and Cuff (use with adaptor)
5082-22	Adult Inflation System with Two-Tube Bag and Cuff (use with adaptor)
5082-23	Large Adult Inflation System with Two-Tube Bag and Cuff (use with adaptor)
5082-45	Thigh, Inflation System with Two-Tube Bag and Cuff (use with adaptor)
1DSL	DuraShock Adapter with one 24" tube and an economy bulb and valve



5085-09
Storage Case, Vinyl



5086-08
Bulb, Standard, Gray, Medium



5087-16
Air Release Valve, Standard



5082-257
Integrated Gauge Adapter

HOME BLOOD PRESSURE MEASUREMENT

BLOOD PRESSURE MANAGEMENT



7052-34

OSZ 5 Digital Blood Pressure System Auto Inflation



7052-40

OSZ 3 Digital Blood Pressure System, Auto Wrist Inflation

OSZ 5 Home Blood Pressure Systems

Accurate, easy-to-use automatic blood pressure systems.

- Easy-to-use one-button operation
- "Fuzzy Logic" smart inflation system is faster and more comfortable
- Large, easy-to-read display shows systolic, diastolic, and pulse measurements simultaneously
- Latex-free

7052-33 OSZ 4 Self Measurement Blood Pressure System, Manual Inflation

7052-34 OSZ 5 Digital Blood Pressure System, Auto Inflation with Adult Cuff, 22-32 cm (8.66 - 12.6 inches)

7052-40 OSZ 3 Digital Blood Pressure System, Auto Wrist Inflation with Universal Cuff (Fits all wrists)

Home Blood Pressure Systems: Accessories (Adult and Large Adult Cuff Only)

7052-35 Adult Cuff for OSZ 5, Digital Blood Pressure System, 22-32cm

7052-36 Large Adult Cuff for OSZ 5, Digital Blood Pressure System, 32 - 42 cm (12.6 - 16.5 in)

7052-38 Voltage Adapter (OSZ 5 only)

7052-38 "L" Connector

7052-39 "T" Connector

FlexiPort® Reusable One-Piece Blood Pressure Cuffs

Long lasting, easy to use, and affordable, Welch Allyn cuffs are the top choice for hospitals and clinics around the world. Why? Because Welch Allyn can offer you virtually every type, size, and configuration of cuffs, along with the ability to simplify inventory and reduce confusion with our innovative FlexiPort connection, which makes every cuff a one- or two-tube cuff.

- Primary sizes are color coded to make it easy to find the proper cuff size
- Easy one-hand attachment/detachment of tubing with FlexiPort makes changing cuffs much faster and easier
- Latex-free material minimizes risk of allergic reactions
- Rotatable port reduces stress to cuff tubing and port, improving patient comfort and cuff durability
- Folded edge reduces risk of cuts and scrapes, maximizing patient comfort
- Meets all latest clinical guidelines for proper fit from AAMI and AHA



FlexiPort Reusable Blood Pressure Cuffs with No Tubes or Connector (Raw Cuff)—One each

REUSE-13	Thigh [size 13]
REUSE-12L	Large Adult Long [size 12L]
REUSE-12	Large Adult [size 12]
REUSE-11L	Adult Long [size 11L]
REUSE-11	Adult [size 11]
REUSE-10	Small Adult [size 10]
REUSE-09	Child [size 9]
REUSE-08	Small Child [size 8]
REUSE-07	Infant [size 7]
REUSE-06	Small Infant [size 6]

FlexiPort Reusable Blood Pressure Cuffs with Screw-Type Connector—One each

REUSE-13-1SC	Thigh [size 13]
REUSE-12L-1SC	Large Adult Long [size 12L]
REUSE-12-1SC	Large Adult [size 12]
REUSE-11L-1SC	Adult Long [size 11L]
REUSE-11-1SC	Adult [size 11]
REUSE-10-1SC	Small Adult [size 10]
REUSE-09-1SC	Child [size 9]



FlexiPort® Adult Reusable Blood Pressure Cuff with No Tubes or Connectors



FlexiPort® Adult Reusable Blood Pressure Cuff with One-Tube Screw-Type Connector

FLEXIPOINT BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT

REUSE-08-1SC Small Child [size 8]

REUSE-07-1SC Infant [size 7]

REUSE-06-1SC Small Infant [size 6]

FlexiPort Reusable Blood Pressure Cuffs with One-Tube Tri-Purpose Connectors—One each

REUSE-13-1TP Thigh [size 13]

REUSE-12L-1TP Large Adult Long [size 12L]

REUSE-12-1TP Large Adult [size 12]

REUSE-11L-1TP Adult Long [size 11L]

REUSE-11-1TP Adult [size 11]

REUSE-10-1TP Small Adult [size 10]

REUSE-09-1TP Child [size 9]

REUSE-08-1TP Small Child [size 8]

REUSE-07-1TP Infant [size 7]

REUSE-06-1TP Small Infant [size 6]

FlexiPort Reusable Blood Pressure Cuffs with One-Tube Locking Type Connectors—One each

REUSE-13-1MQ Thigh [size 13]

REUSE-12L-1MQ Large Adult Long [size 12L]

REUSE-12-1MQ Large Adult [size 12]

REUSE-11L-1MQ Adult Long [size 11L]

REUSE-11-1MQ Adult [size 11]

REUSE-10-1MQ Small Adult [size 10]

REUSE-09-1MQ Child [size 9]

REUSE-08-1MQ Small Child [size 8]

REUSE-07-1MQ Infant [size 7]

REUSE-06-1MQ Small Infant [size 6]

FlexiPort Reusable Blood Pressure Cuffs with One-Tube Bayonet Connectors—One each

REUSE-13-1HP Thigh [size 13]

REUSE-12L-1HP Large Adult Long [size 12L]

REUSE-12-1HP Large Adult [size 12]

REUSE-11L-1HP Adult Long [size 11L]

REUSE-11-1HP	Adult [size 11]
REUSE-10-1HP	Small Adult [size 10]
REUSE-09-1HP	Child [size 9]
REUSE-08-1HP	Small Child [size 8]
REUSE-07-1HP	Infant [size 7]
REUSE-06-1HP	Small Infant [size 6]

FlexiPort Reusable Blood Pressure Cuffs with Two-Tube Screw-type Connectors—One each

REUSE-13-2SC	Thigh [size 13]
REUSE-12L-2SC	Large Adult Long [size 12L]
REUSE-12-2SC	Large Adult [size 12]
REUSE-11L-2SC	Adult Long [size 11L]
REUSE-11-2SC	Adult [size 11]
REUSE-10-2SC	Small Adult [size 10]
REUSE-09-2SC	Child [size 9]
REUSE-08-2SC	Small Child [size 8]
REUSE-07-2SC	Infant [size 7]

FlexiPort Reusable Blood Pressure Cuffs with Two-Tube Locking-type Connectors—One each

REUSE-13-2MQ	Thigh [size 13]
REUSE-12L-2MQ	Large Adult Long [size 12L]
REUSE-12-2MQ	Large Adult [size 12]
REUSE-11L-2MQ	Adult Long [size 11L]
REUSE-11-2MQ	Adult [size 11]
REUSE-10-2MQ	Small Adult [size 10]
REUSE-09-2MQ	Child [size 9]
REUSE-08-2MQ	Small Child [size 8]
REUSE-07-2MQ	Infant [size 7]



FlexiPort® Adult Reusable Blood Pressure Cuff with Two-Tube Screw-Type Connector

FLEXIPOINT BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT

FlexiPort Reusable Blood Pressure Cuffs with Two-Tube Tri-Purpose Connector & Empty Tube—One each

REUSE-13-2TP	Thigh [size 13]
REUSE-12L-2TP	Large Adult Long [size 12L]
REUSE-12-2TP	Large Adult [size 12]
REUSE-11L-2TP	Adult Long [size 11L]
REUSE-11-2TP	Adult [size 11]
REUSE-10-2TP	Small Adult [size 10]
REUSE-09-2TP	Child [size 9]
REUSE-08-2TP	Small Child [size 8]
REUSE-07-2TP	Infant [size 7]

FlexiPort Reusable Blood Pressure Cuffs with Inflation System (bulb and valve)

REUSE-13-2BV	Thigh [size 13]
REUSE-12L-2BV	Large Adult Long [size 12L]
REUSE-12-2BV	Large Adult [size 12]
REUSE-11L-2BV	Adult Long [size 11L]
REUSE-11-2BV	Adult [size 11]
REUSE-10-2BV	Small Adult [size 10]
REUSE-09-2BV	Child [size 9]
REUSE-08-2BV	Small Child [size 8]
REUSE-07-2BV	Infant [size 7]

FlexiPort® Soft Disposable One-Piece Blood Pressure Cuffs

FlexiPort disposable blood pressure cuffs are the top choice for hospitals and clinics around the world. That's because we offer virtually every type of cuff an institution may need, along with the ability to simplify inventory like no other cuff provider with our patented FlexiPort connection.

- Every cuff can be used as a one- or two-tube cuff thanks to patented FlexiPort technology
- Every cuff size, from small infant through thigh cuffs
- Choose from either soft or vinyl fabric to fit anyone's preference
- Color-coded for fast, accurate selection of the correct cuff size
- Rotatable port reduces stress to cuff tubing and port, improving patient comfort and cuff durability
- Meets all the latest clinical guidelines for proper fit from AAMI and AHA



FlexiPort Soft Disposable Blood Pressure Cuffs, No Tubes or Connectors (Raw Cuff) (20/case)

SOFT-13 Thigh [size 13]

SOFT-12L Large Adult Long [size 12L]

SOFT-12 Large Adult [size 12]

SOFT-11L Adult Long [size 11L]

SOFT-11 Adult [size 11]

SOFT-10 Small Adult [size 10]

SOFT-09 Child [size 9]

SOFT-08 Small Child [size 8]

SOFT-07 Infant [size 7]

SOFT-06 Small Infant [size 6]



FlexiPort® Soft Disposable Blood Pressure Cuff (no tubes or connectors)

FlexiPort Soft Disposable Blood Pressure Cuffs, One-Tube Screw-Type Connector (20/case)

SOFT-13-1SC Thigh [size 13]

SOFT-12L-1SC Large Adult Long [size 12L]

SOFT-12-1SC Large Adult [size 12]

SOFT-11L-1SC Adult Long [size 11L]

SOFT-11-1SC Adult [size 11]

SOFT-10-1SC Small Adult [size 10]

SOFT-09-1SC Child [size 9]

SOFT-08-1SC Small Child [size 8]



FlexiPort® Soft Disposable Blood Pressure Cuffs, One-Tube Screw-Type Connector

FLEXIPOINT BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT

SOFT-07-1SC Infant [size 7]

SOFT-06-1SC Small Infant [size 6]

FlexiPort Soft Disposable Blood Pressure Cuffs, One-Tube Tri-Purpose Connector (20/case)

SOFT-13-1TP Thigh [size 13]

SOFT-12L-1TP Large Adult Long [size 12L]

SOFT-12-1TP Large Adult [size 12]

SOFT-11L-1TP Adult Long [size 11L]

SOFT-11-1TP Adult [size 11]

SOFT-10-1TP Small Adult [size 10]

SOFT-09-1TP Child [size 9]

SOFT-08-1TP Small Child [size 8]

SOFT-07-1TP Infant [size 7]

SOFT-06-1TP Small Infant [size 6]

FlexiPort Soft Disposable Blood Pressure Cuffs, One-Tube Locking-Type Connector (20/case)

SOFT-13-1MQ Thigh [size 13]

SOFT-12L-1MQ Large Adult Long [size 12L]

SOFT-12-1MQ Large Adult [size 12]

SOFT-11L-1MQ Adult Long [size 11L]

SOFT-11-1MQ Adult [size 11]

SOFT-10-1MQ Small Adult [size 10]

SOFT-09-1MQ Child [size 9]

SOFT-08-1MQ Small Child [size 8]

SOFT-07-1MQ Infant [size 7]

SOFT-06-1MQ Small Infant [size 6]

FlexiPort Soft Disposable Blood Pressure Cuffs, One-Tube Bayonet-Type Connector (20/case)

SOFT-13-1HP Thigh [size 13]

SOFT-12L-1HP Large Adult Long [size 12L]

SOFT-12-1HP Large Adult [size 12]

SOFT-11L-1HP Adult Long [size 11L]

SOFT-11-1HP Adult [size 11]

SOFT-10-1HP Small Adult [size 10]

SOFT-09-1HP	Child [size 09]
SOFT-08-1HP	Small Child [size 8]
SOFT-07-1HP	Infant [size 7]
SOFT-06-1HP	Small Infant [size 6]

FlexiPort Soft Disposable Blood Pressure Cuffs, Two-Tube Screw Type Connectors (20/case)

SOFT-13-2SC	Thigh [size 13]
SOFT-12L-2SC	Large Adult Long [size 12L]
SOFT-12-2SC	Large Adult [size 12]
SOFT-11L-2SC	Adult Long [size 11L]
SOFT-11-2SC	Adult [size 11]
SOFT-10-2SC	Small Adult [size 10]
SOFT-09-2SC	Child [size 9]
SOFT-08-2SC	Small Child [size 8]
SOFT-07-2SC	Infant [size 7]

FlexiPort Soft Disposable Blood Pressure Cuffs with Two-Tube Locking Type Connector (20/case)

SOFT-13-2MQ	Thigh [size 13]
SOFT-12L-2MQ	Large Adult Long [size 12L]
SOFT-12-2MQ	Large Adult [size 12]
SOFT-11L-2MQ	Adult Long [size 11L]
SOFT-11-2MQ	Adult [size 11]
SOFT-10-2MQ	Small Adult [size 10]
SOFT-09-2MQ	Child [size 09]
SOFT-08-2MQ	Small Child [size 08]
SOFT-07-2MQ	Infant [size 7]

FlexiPort Soft Disposable Blood Pressure Cuffs with Two-Tube Tri-Purpose Connector/Empty Tube (20/case)

SOFT-13-2TP	Thigh [size 13]
SOFT-12L-2TP	Large Adult Long [size 12L]
SOFT-12-2TP	Large Adult [size 12]
SOFT-11L-2TP	Adult Long [size 11L]
SOFT-11-2TP	Adult [size 11]



FlexiPort Soft Disposable Blood Pressure Cuffs with Two-Tube Screw-Type Connector

FLEXIPOINT BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT

SOFT-10-2TP Small Adult [size 10]

SOFT-09-2TP Child [size 9]

SOFT-08-2TP Small Child [size 8]

SOFT-07-2TP Infant [size 7]

FlexiPort Soft Disposable Blood Pressure Cuffs with Two-Tube Inflation System (20/case)

SOFT-13-2BV Thigh [size 13]

SOFT-12L-2BV Large Adult Long [size 12L]

SOFT-12-2BV Large Adult [size 12]

SOFT-11L-2BV Adult Long [size 11L]

SOFT-11-2BV Adult [size 11]

SOFT-10-2BV Small Adult [size 10]

SOFT-09-2BV Child [size 09]

SOFT-08-2BV Small Child [size 08]

SOFT-07-2BV Infant [size 7]

FlexiPort® Vinyl Disposable One-Piece Blood Pressure Cuffs

Welch Allyn FlexiPort disposable blood pressure cuffs are the top choice for hospitals and clinics around the world. That's because we offer virtually every type of cuff an institution may need, along with the ability to simplify inventory like no other cuff provider with our patented FlexiPort connection.

- Every cuff can be used as a one- or two-tube cuff thanks to patented FlexiPort technology
- Every cuff size, from small infant through thigh cuffs
- Choose from either soft or vinyl fabric to fit anyone's preference
- Color-coded for fast, accurate selection of the correct cuff size
- Rotatable port reduces stress to cuff tubing and port, improving patient comfort and cuff durability
- Meets all the latest clinical guidelines for proper fit from AAMI and AHA



FlexiPort Vinyl Blood Pressure Cuffs, No Tubes or Connectors (Raw Cuff) (20/case)

VINYL-13	Thigh [size 13]
VINYL-12L	Large Adult Long [size 12L]
VINYL-12	Large Adult [size 12]
VINYL-11L	Adult Long [size 11L]
VINYL-11	Adult [size 11]
VINYL-10	Small Adult [size 10]
VINYL-09	Child [size 9]
VINYL-08	Small Child [size 8]
VINYL-07	Infant [size 7]
VINYL-06	Small Infant [size 6]



FlexiPort® Vinyl Disposable Blood Pressure Cuff (no tubes or connectors)

FlexiPort Vinyl Blood Pressure Cuffs, One-Tube Screw-Type Connector (20/case)

VINYL-13-1SC	Thigh [size 13]
VINYL-12L-1SC	Large Adult Long [size 12L]
VINYL-12-1SC	Large Adult [size 12]
VINYL-11L-1SC	Adult Long [size 11L]
VINYL-11-1SC	Adult [size 11]
VINYL-10-1SC	Small Adult [size 10]
VINYL-09-1SC	Child [size 9]
VINYL-08-1SC	Small Child [size 8]



FlexiPort® Vinyl Disposable Blood Pressure Cuffs, One-Tube Screw-Type Connector

FLEXIPOINT BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT

VINYL-07-1SC	Infant [size 7]
VINYL-06-1SC	Small Infant [size 6]

FlexiPort Vinyl Disposable Blood Pressure Cuffs, One-Tube Locking-Type Connector (20/case)

VINYL-13-1MQ	Thigh [size 13]
VINYL-12L-1MQ	Large Adult Long [size 12L]
VINYL-12-1MQ	Large Adult [size 12]
VINYL-11L-1MQ	Adult Long [size 11L]
VINYL-11-1MQ	Adult [size 11]
VINYL-10-1MQ	Small Adult [size 10]
VINYL-09-1MQ	Child [size 9]
VINYL-08-1MQ	Small Child [size 8]
VINYL-07-1MQ	Infant [size 7]
VINYL-06-1MQ	Small Infant [size 6]

FlexiPort Vinyl Disposable Blood Pressure Cuffs, One-Tube Bayonet-Type Connector (20/case)

VINYL-13-1HP	Thigh [size 13]
VINYL-12L-1HP	Large Adult Long [size 12L]
VINYL-12-1HP	Large Adult [size 12]
VINYL-11L-1HP	Adult Long [size 11L]
VINYL-11-1HP	Adult [size 11]
VINYL-10-1HP	Small Adult [size 10]
VINYL-09-1HP	Child [size 09]
VINYL-08-1HP	Small Child [size 8]
VINYL-07-1HP	Infant [size 7]
VINYL-06-1HP	Small Infant [size 6]

FlexiPort Vinyl Disposable Blood Pressure Cuffs, Two-Tube Screw Type Connectors (20/case)

VINYL-13-2SC	Thigh [size 13]
VINYL-12L-2SC	Large Adult Long [size 12L]
VINYL-12-2SC	Large Adult [size 12]
VINYL-11L-2SC	Adult Long [size 11L]
VINYL-11-2SC	Adult [size 11]
VINYL-10-2SC	Small Adult [size 10]
VINYL-09-2SC	Child [size 9]
VINYL-08-2SC	Small Child [size 8]
VINYL-07-2SC	Infant [size 7]

FlexiPort Vinyl Disposable Blood Pressure Cuffs with Two-Tube Locking-Type Connector (20/case)

VINYL-13-2MQ	Thigh [size 13]
VINYL-12L-2MQ	Large Adult Long [size 12L]
VINYL-12-2MQ	Large Adult [size 12]
VINYL-11L-2MQ	Adult Long [size 11L]
VINYL-11-2MQ	Adult [size 11]
VINYL-10-2MQ	Small Adult [size 10]
VINYL-09-2MQ	Child [size 09]
VINYL-08-2MQ	Small Child [size 08]
VINYL-07-2MQ	Infant [size 7]

FlexiPort Vinyl Disposable Blood Pressure Cuffs with Two-Tube Inflation System (20/case)

VINYL-13-2BV	Thigh [size 13]
VINYL-12L-2BV	Large Adult Long [size 12L]
VINYL-12-2BV	Large Adult [size 12]
VINYL-11L-2BV	Adult Long [size 11L]
VINYL-11-2BV	Adult [size 11]
VINYL-10-2BV	Small Adult [size 10]
VINYL-09-2BV	Child [size 09]
VINYL-08-2BV	Small Child [size 08]
VINYL-07-2BV	Infant [size 7]



FlexiPort Vinyl Disposable Blood Pressure Cuffs with Two-tube Screw Type Connector

FLEXIPOINT ECOCUFF BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT



FlexiPort® EcoCuff™ Single-Patient-Use Blood Pressure Cuffs

Featuring Welch Allyn FlexiPort technology and a number of clinical enhancements that facilitate proper sizing, the EcoCuff helps you maximize value while minimizing cross-contamination related to BP cuffs. And because it's capable of being recycled and built using fewer, lighter materials than traditional disposable blood pressure cuffs, the EcoCuff can reduce the volume of waste in your facility when used in a one-cuff-per-patient model.

- FlexiPort® single-point connection standardization system lets you use one cuff in virtually every room and on virtually every device you have in house
- Lifecycle analysis demonstrates over 50% less environmental impact than other FlexiPort disposable cuffs
- Durable enough to last for a complete patient stay
- Unique slot fitting system helps ensure patient is receiving the proper-sized cuff
- Four sizes fit over 90% of the U.S. patient population over one year old, making it easy to find the proper cuff and fewer part numbers simplify inventory management
- Meets AAMI standards

FlexiPort EcoCuff Single-Patient-Use Blood Pressure Cuff (Raw Cuff) (20 cuffs per bag, 5 bags per case)

Eco-09 Child, 15-21 cm

Eco-10 Small Adult, 20-28 cm

Eco-11 Adult, 27-38 cm

Eco-12 Large Adult, 33-45 cm



The EcoCuff is capable of being recycled. Although the practice is not yet widespread, a number of hospitals have implemented the practice of recycling medical-use polypropylene plastics.

FlexiPort® Accessories

1-BARB One-Tube FlexiPort Double Hose Barb Connector (8")

1-SC One-Tube FlexiPort Screw Connector (8")

1-MQ One-Tube FlexiPort Locking Connector (8")

1-HP One-Tube FlexiPort Bayonet Connector (8")

1-TPS One-Tube FlexiPort Tri-Purpose Connector
(24" for hand aneroid devices)

1-TPSE One-Tube FlexiPort Tri-Purpose Connector
(24" for European hand aneroid devices)

1-TPL One-Tube FlexiPort Tri-Purpose Connector
(31" for hand aneroid devices)

1-TPLE One-Tube FlexiPort Tri-Purpose Connector
(31" for European hand aneroid devices)

2-BARB Two-Tube FlexiPort Double Hose Barb Connector (8")

2-SC Two-Tube FlexiPort Screw Connector (8")

2-SCC Two-Tube FlexiPort Screw Connector (8") with Clamp on Tube

2-MQ Two-Tube FlexiPort Locking Connector (8")

2-MF Two-Tube FlexiPort Mated Locking Connector (8")

2-TPLE Two-Tube FlexiPort Tri-Purpose Connector
(31" for European wall, mobile and pocket aneroids)

2-TPSE Two-Tube FlexiPort Tri-Purpose Connector
(24" for European wall, mobile and pocket aneroids)

2-TP Two-Tube FlexiPort Tri-Purpose Connector
(8" and 13" for wall, mobile and pocket aneroids)

2-SRC Two-Tube FlexiPort Tri-Purpose Connector
(8" and 13" for wall, mobile and pocket aneroids)



1-BARB



1-SC



1-MQ



1-HP



1-TPS 1-TPSE 1-TPL 1-TPLE



2-BARB



2-SC



2-MQ



2-MF



2-TPLE, 2-TPSE



2-TP, 2-SRC

FLEXIPOINT BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT



2-BVL



2-BVD

5082-188



5082-189



PORT-1



PORT-2



5082-159

Wash Cap that Tethers to FlexiPort on Cuffs



1SRCL, 1SRCS



1DSAL, 1DSAS



1DSLCL

2-BVL Two-Tube FlexiPort Tri-Purpose Connector (8" and 13") with Large Bulb

2-BVLC Two-Tube FlexiPort Tri-Purpose Connector (8" and 13") with Large Bulb and Clamp

2-BVS Two-Tube FlexiPort Tri-Purpose Connector (8" and 13") with Small Bulb

2-BVSC Two-Tube FlexiPort Tri-Purpose Connector (8" and 13") with Small Bulb and Clamp

2-BVD Two-Tube FlexiPort Tri-Purpose Connector (8" and 13") with Econo Bulb

2-BVDC Two-Tube FlexiPort Tri-Purpose Connector (8" and 13") with Economy Bulb and Clamp

5082-188 8" FlexiPort Tubing

5082-189 13" FlexiPort Tubing

PORT-1 One-Tube FlexiPort Fitting without any Tubes

PORT-2 Two-Tube FlexiPort Fitting without any Tubes

5082-159 Wash Cap that Tethers to FlexiPort on Cuffs

1SRCL DuraShock Classic Adapter with 24" Tube and Large Bulb with Premium Valve

1SRCS DuraShock Classic Adapter with 24" Tube and Small Bulb with Premium Valve

1DSAL DuraShock Adapter with 24" Tube and Large Bulb with Premium Valve

1DSAS DuraShock Adapter with 24" Tube and Small Bulb with Premium Valve

1DSLCL DuraShock Adapter with 24" Tube and Economy Bulb and Valve

Neonatal Soft Disposable One-Piece Blood Pressure Cuffs with Luer Slip Connectors (10 pack)

Reliable, easy-to-use, and affordable.

- Soft, comfortable DuPont® Sontara® spunlaced fabric is durable and 100% latex-free
- Range/index markings on outside of cuff make sizing easy
- Color coded for easy size identification
- Flat port tubing connectors help keep tubes out of the way
- Connectors to fit virtually any brand of NIBP monitors



One-Tube Neonatal Soft Disposable One-Piece Blood Pressure Cuffs with Luer Slip Connectors

5082-101-1 Neonate [size 1] 1 Tube, 10 Pack

5082-102-1 Neonate [size 2] 1 Tube, 10 Pack

5082-103-1 Neonate [size 3] 1 Tube, 10 Pack

5082-104-1 Neonate [size 4] 1 Tube, 10 Pack

5082-105-1 Neonate [size 5] 1 Tube, 10 Pack

Two-Tube Neonatal Soft Disposable One-Piece Blood Pressure Cuffs with Luer Slip Connectors

5082-101-2 Neonate [size 1] 2 Tube, 10 Pack

5082-102-2 Neonate [size 2] 2 Tube, 10 Pack

5082-103-2 Neonate [size 3] 2 Tube, 10 Pack

5082-104-2 Neonate [size 4] 2 Tube, 10 Pack

5082-105-2 Neonate [size 5] 2 Tube, 10 Pack

TRIMLINE BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT



Trimline™ Reusable Blood Pressure Cuffs

A traditional style blood pressure cuff with enhanced design to facilitate accurate readings.

- Color coding helps clinicians quickly identify the proper size cuff
- Completely latex-free with antimicrobial agents applied to help prevent bacterial growth and mildew on the cuff fabric
- Connects to all manual and NIBP monitors
- Exterior Pre-Gaged® range guide and BALANCED® bladder design help to ensure accurate readings by facilitating proper arterial compression

Reusable One-Tube Cuffs with Female Slip Luer Connector (5 per box / Each)



3600S/1900S	Thigh
3602XLS/1902XLS	Large Adult Long
3602S/1902S	Large Adult
3603XLS/1903XLS	Adult Long
3603S/1903S	Adult
3604S/1904S	Small Adult
3605S/1905S	Child
3606S/1906S	Infant

Reusable One-Tube Cuffs with Female Subminiature Connector (5 per box / Each)

3600SMQ/1900SMQ	Thigh
3602XLSMQ/1902XLSMQ	Large Adult Long
3602SMQ/1902SMQ	Large Adult
3603XLSMQ/1903XLSMQ	Adult Long
3603SMQ/1903SMQ	Adult
3604SMQ/1904SMQ	Small Adult
3605SMQ/1905SMQ	Child
3606SMQ/1906SMQ	Infant

Reusable One-Tube Cuffs with Male Bayonet/Rectus Connector (5 per box / Each)

3600SHP/1900SHP	Thigh
3602XLSHP/1902XLSHP	Large Adult Long
3602SHP/1902SHP	Large Adult
3603XLSHP/1903XLSHP	Adult Long
3603SHP/1903SHP	Adult
3604SHP/1904SHP	Small Adult
3605SHP/1905SHP	Child
3606SHP/1906SHP	Infant

Reusable One-Tube Cuffs with Male Subminiature Connector (5 per box / Each)

36700	Thigh
36709	Large Adult Long
36702	Large Adult
36708	Adult Long
36703	Adult
36704	Small Adult
36705	Child
36706	Infant

Reusable Two-Tube Cuffs with one Female Slip Luer Connector and Bare Tube (5 per box / Each)

3600/1900	Thigh
3602XL/1902XL	Large Adult Long
3602/1902	Large Adult
3603XL/1903XL	Adult Long
3603/1903	Adult
3604/1904	Small Adult
3605/1905	Child
3606/1906	Infant

TRIMLINE BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT



Reusable Two-Tube Cuffs with Female Slip Luer Connectors (5 per box / Each)

3600F/1900F	Thigh
3602XLF/1902XLF	Large Adult Long
3602F/1902F	Large Adult
3603XLF/1903XLF	Adult Long
3603F/1903F	Adult
3604F/1904F	Small Adult
3605F/1905F	Child
3606F/1906F	Infant

Reusable Two-Tube Cuffs with Screw Connectors (5 per box / Each)

3600CK/1900CK	Thigh
3602XLCK/1902XLCK	Large Adult Long
3602CK/1902CK	Large Adult
3603XLCK/1903XLCK	Adult Long
3603CK/1903CK	Adult
3604CK/1904CK	Small Adult
3605CK/1905CK	Child
3606CK/1906CK	Infant

Reusable Two-Tube Cuffs with Female Subminiature Connectors (5 per box / Each)

3600MQ/1900MQ	Thigh
3602XLMQ/1902XLMQ	Large Adult Long
3602MQ/1902MQ	Large Adult
3603XLMQ/1903XLMQ	Adult Long
3603MQ/1903MQ	Adult
3604MQ/1904MQ	Small Adult
3605MQ/1905MQ	Child
3606MQ/1906MQ	Infant

Reusable Two-Tube Cuffs with Mated Subminiature Connectors (5 per box / Each)

3600DX/1900DX	Thigh
3602XLDX/1902XLDX	Large Adult Long
3602DX/1902DX	Large Adult
3603XLDX/1903XLDX	Adult Long
3603DX/1903DX	Adult
3604DX/1904DX	Small Adult
3605DX/1905DX	Child
3606DX/1906DX	Infant

Reusable Two-Tube Cuffs with Bulb and Valve and Female Slip Luer connector (5 per box / Each)

3300/1600	Thigh
3302XL/1602XL	Large Adult Long
3302/1602	Large Adult
3303XL/1603XL	Adult Long
3303/1603	Adult
3304/1604	Small Adult
3305/1605	Child
3306/1606	Infant

TRIMLINE BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT



Trimline™ Soft Disposable Blood Pressure Cuffs

For facilities looking for a traditional blood pressure cuff design, Trimline disposable cuffs offer a comfortable and economical solution.

- Soft against the skin for patient comfort
- Available in 10 sizes with several connector configurations, so you can standardize on our cuffs to work on virtually any blood pressure device
- Color-coded printing on the fabric allows easy size identification
- Long cuffs are approximately 4-5" longer than standard cuffs to accommodate very large arms or bariatric patients
- Exterior Pre-Gaged® range guide and BALANCED® bladder design help to ensure accurate readings by facilitating proper arterial compression

Soft Disposable One-Tube Cuffs with Bare Tube

39001	Thigh
39240	Large Adult Long
39124	Large Adult
39042	Adult Long
39004	Adult
39027	Small Adult
39083	Child
39055	Child Long
39117	Small Child
39239	Infant

Soft Disposable One-Tube Cuffs with Screw Connector

39076	Thigh
39075	Large Adult Long
39074	Large Adult
39093	Adult Long
39073	Adult
39092	Small Adult
39091	Child
39069	Child Long
39135	Small Child
39090	Infant



39004
One-tube Adult

Soft Disposable One-Tube Cuffs with Male Bayonet/Rectus Connector

39015	Thigh
39060	Large Adult Long
39230	Large Adult
39019	Adult Long
39048	Adult
39085	Small Adult
39100	Child
39173	Child Long
39024	Small Child
39045	Infant

Soft Disposable One-Tube Cuffs with Female Subminiature Connector

39236	Thigh
39237	Large Adult Long
39039	Large Adult
39248	Adult Long
39080	Adult
39238	Small Adult
39007	Child
39108	Child Long
39058	Small Child
39026	Infant

Soft Disposable One-Tube Cuffs with Male Subminiature Connector

39102	Thigh
39158	Large Adult Long
39261	Large Adult
39270	Adult Long
39105	Adult
39242	Small Adult
39044	Child
39145	Child Long
39020	Small Child
39127	Infant

Soft Disposable Two-Tube Cuffs with Bare Tubes

39148	Thigh
39047	Large Adult Long
39103	Large Adult
39028	Adult Long
39246	Adult
39107	Small Adult
39154	Child
39070	Child Long
39165	Small Child
39008	Infant

Soft Disposable Two-Tube Cuffs with Screw Connectors

39072	Thigh
39071	Large Adult Long
39014	Large Adult
39053	Adult Long
39023	Adult
39012	Small Adult
39011	Child
39038	Child Long
39128	Small Child
39010	Infant

Soft Disposable Two-Tube Cuffs with Screw Connectors and Clamp

39218	Thigh
39189	Large Adult Long
39066	Large Adult
39065	Adult Long
39064	Adult
39063	Small Adult
39062	Child
39129	Child Long
39077	Small Child
39061	Infant

Soft Disposable Two-Tube Cuffs with Female Subminiature Connectors

39156	Thigh
39018	Large Adult Long
39149	Large Adult
39081	Adult Long
39268	Adult
39213	Small Adult
39249	Child
39134	Child Long
39046	Small Child
39163	Infant

Soft Disposable Two-Tube Cuffs with Female Slip Luer Connectors

39164	Thigh
39003	Large Adult Long
39157	Large Adult
39267	Adult Long
39041	Adult
39169	Small Adult
39250	Child
39214	Child Long
39170	Small Child
39052	Infant

Soft Disposable Two-Tube Cuffs with Mated Subminiature Connectors

39037	Thigh
39036	Large Adult Long
39035	Large Adult
39034	Adult Long
39033	Adult
39032	Small Adult
39031	Child
39168	Child Long
39151	Small Child
39030	Infant

TRIMLINE BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT

Soft Disposable Two-Tube Cuffs with Bulb, Valve, Clamp and Bare tube

39029	Thigh
39125	Large Adult Long
39254	Large Adult
39126	Adult Long
39167	Adult
39050	Small Adult
39115	Child
39200	Child Long
39188	Small Child
39171	Infant

Soft Disposable Two-Tube Cuffs with Bulb, Valve, Clamp and One Screw Connector

39131	Thigh
39229	Large Adult Long
39021	Large Adult
39172	Adult Long
39122	Adult
39121	Small Adult
39120	Child
39257	Child Long
39235	Small Child
39119	Infant

Soft Disposable Two-Tube Cuffs with Bulb, Valve, Clamp and One Bayonet/Rectus Connector

39054	Thigh
39094	Large Adult Long
39272	Large Adult
39160	Adult Long
39199	Adult
39049	Small Adult
39234	Child
39243	Child Long
39216	Small Child
39155	Infant

Trimline™ Vinyl Disposable Blood Pressure Cuffs

Wipeable vinyl cuffs are optimal for use in the Emergency Room and other areas where fluids present cleaning concerns.

- Constructed of moisture-resistant, wipeable material
- Eight sizes of cuffs available in configurations most often used in ED and other patient care areas
- Color-coded, alcohol resistant ink allows quick identification by size
- Available with a bulb and valve or a pocket aneroid gauge for manual blood pressure monitoring situations
- Two sizes with a longer length are available to accommodate very large arms or bariatric patients
- Exterior Pre-Gaged® range guide and BALANCED® bladder design help to ensure accurate readings by facilitating proper arterial compression



99859
One-Tube Adult

Vinyl Disposable One-Tube Cuffs with Screw Connector

99865	Thigh
99862	Large Adult Long
99863	Large Adult
99858	Adult Long
99859	Adult
99855	Small Adult
99853	Child
99851	Infant

Vinyl Disposable One-Tube Cuffs with Male Bayonet/Rectus Connector

99874	Thigh
99872	Large Adult Long
99871	Large Adult
99870	Adult Long
99869	Adult
99868	Small Adult
99867	Child
99866	Infant

TRIMLINE BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT



99825
Two-Tube Adult

Vinyl Disposable Two-Tube Cuffs with Bare Tubes

99833	Thigh
99804	Large Adult Long
99837	Large Adult
99811	Adult Long
99825	Adult
99820	Small Adult
99816	Child
99842	Infant

Vinyl Disposable Two-Tube Cuffs with Screw Connectors

99864	Thigh
99861	Large Adult Long
99860	Large Adult
99857	Adult Long
99856	Adult
99854	Small Adult
99852	Child
99850	Infant

Vinyl Disposable Two-Tube Cuffs with Screw Connectors and Clamp

99884	Thigh
99883	Large Adult Long
99881	Large Adult
99880	Adult Long
99879	Adult
99878	Small Adult
99877	Child
99876	Infant

Vinyl Disposable Two-tube Cuffs with Bulb, Valve, Clamp and One Screw Connector

99885	Large Adult
99886	Adult Long
99882	Adult

Trimline™ Disposable Isolation Blood Pressure Cuffs

Reduce the risk of cross-contamination with a cuff designed to be used on a single patient in isolation situations.

- Bright yellow fabric for fast “Isolation Use” identification
- Soft against the skin for patient comfort
- Available in 8 sizes with several connector configurations
- Exterior Pre-Gaged® range guide and BALANCED® bladder design help to ensure accurate readings by facilitating proper arterial compression



39203
Adult Two-Tube

Isolation One-Tube Cuffs with Male Bayonet/Rectus Connector

39221	Thigh
39227	Large Adult Long
39222	Large Adult
39228	Adult Long
39223	Adult
39224	Small Adult
39225	Child
39226	Infant

Isolation Two-Tube Cuffs with Screw Connectors with Clamp

39201	Thigh
39207	Large Adult Long
39202	Large Adult
39208	Adult Long
39203	Adult
39204	Small Adult
39205	Child
39206	Infant

Two-Tube Cuffs with Bulb, Valve, Clamp and Bare Tube

39212	Large Adult
39211	Adult Long
39209	Adult

TRIMLINE BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT



Trimline Soft Neonatal One-Piece Blood Pressure Cuffs

Accuracy and comfort are crucial for your tiniest patients. Welch Allyn offers a choice of neonatal blood pressure cuffs in soft fabric or supple vinyl material.

- Cuffs feature playful graphics and are numbered 1-5 for easy size identification
- The BALANCED® bladder design helps to ensure accurate readings by facilitating proper arterial compression on both right and left arms

One-Tube Trimline Soft Neonatal One-Piece Blood Pressure Cuffs with Luer Slip Connectors

- 39190** Neonate [size 1] 1 Tube, 40 per Case
- 39191** Neonate [size 2] 1 Tube, 40 per Case
- 39192** Neonate [size 3] 1 Tube, 40 per Case
- 39193** Neonate [size 4] 1 Tube, 40 per Case
- 39194** Neonate [size 5] 1 Tube, 40 per Case

Two-Tube Trimline Soft Neonatal One-Piece Blood Pressure Cuffs with Luer Slip Connectors

- 39110** Neonate [size 1] 2 Tube, 40 per Case
- 39111** Neonate [size 2] 2 Tube, 40 per Case
- 39112** Neonate [size 3] 2 Tube, 40 per Case
- 39113** Neonate [size 4] 2 Tube, 40 per Case
- 39115** Neonate [size 5] 2 Tube, 40 per Case

Two-Tube Trimline Soft Neonatal One-Piece Blood Pressure Cuffs with Luer Slip Connectors and Clamp

- 39140** Neonate [size 1] 2 Tube, 40 per Case
- 39141** Neonate [size 2] 2 Tube, 40 per Case
- 39142** Neonate [size 3] 2 Tube, 40 per Case
- 39143** Neonate [size 4] 2 Tube, 40 per Case
- 39144** Neonate [size 5] 2 Tube, 40 per Case

Trimline Vinyl Neonatal One-Piece Blood Pressure Cuffs with Luer Slip Connectors

One-Tube Trimline Vinyl Neonatal One-Piece Blood Pressure Cuffs with Luer Slip Connectors

99801	Neonate [size 1]	1 Tube, 40 per Case
99890	Neonate [size 2]	1 Tube, 40 per Case
99729	Neonate [size 3]	1 Tube, 40 per Case
99848	Neonate [size 4]	1 Tube, 40 per Case
99750	Neonate [size 5]	1 Tube, 40 per Case

One-Tube Trimline Vinyl Neonatal One-Piece Blood Pressure Cuffs with Micro Luer Slip Connectors

99789	Neonate [size 1]	1 Tube, 40 per Case
99884	Neonate [size 2]	1 Tube, 40 per Case
99741	Neonate [size 3]	1 Tube, 40 per Case
99836	Neonate [size 4]	1 Tube, 40 per Case
99707	Neonate [size 5]	1 Tube, 40 per Case

Two-Tube Trimline Vinyl Neonatal One-Piece Blood Pressure Cuffs with Luer Slip Connectors

99760	Neonate [size 1]	2 Tube, 40 per Case
99818	Neonate [size 2]	2 Tube, 40 per Case
99873	Neonate [size 3]	2 Tube, 40 per Case
99875	Neonate [size 4]	2 Tube, 40 per Case
99722	Neonate [size 5]	2 Tube, 40 per Case



TRIMLINE BLOOD PRESSURE CUFFS

BLOOD PRESSURE MANAGEMENT

Two-Tube Trimline Vinyl Neonatal One-Piece Blood Pressure Cuffs with Micro Luer Slip Connectors

99832 Neonate [size 1] 2 Tube, 40 per Case

99721 Neonate [size 2] 2 Tube, 40 per Case

99806 Neonate [size 3] 2 Tube, 40 per Case

99714 Neonate [size 4] 2 Tube, 40 per Case

99824 Neonate [size 5] 2 Tube, 40 per Case

One-Tube Blood Pressure Cuffs and Neoprene Bladders

- For use with hand aneroids
- Heavyweight nylon cuff
- Latex-free

5082-07 Newborn Blood Pressure Cuff and One-Tube Bladder

5082-03 Infant Blood Pressure Cuff and One-Tube Bladder

5082-42 Child Blood Pressure Cuff and One-Tube Bladder

5082-62 Child-Print Blood Pressure Cuff and One-Tube Bladder

5082-43 Adult Blood Pressure Cuff and One-Tube Bladder

5082-44 Large Adult Blood Pressure Cuff and One-Tube Bladder

5082-77 Thigh Blood Pressure Cuff and One-Tube Bladder



5082-43
Adult Blood Pressure Cuff and
One-Tube Bladder



5082-200
Blood Pressure Connector
(Tri-Purpose Type)

One-Tube Neoprene Bladders

- Latex-free neoprene for enhanced safety

5089-37 Black Neoprene Newborn Inflation Bladder

5089-25 Black Neoprene Infant Inflation Bladder

5089-04 Black Neoprene Child Inflation Bladder

5089-03 Black Neoprene Adult Inflation Bladder

5089-21 Black Neoprene Large Adult Inflation Bladder

5089-18 Black Neoprene Thigh Inflation Bladder



5089-03
Black Neoprene Adult
Inflation Bladder

TWO-PIECE CUFFS AND BLADDERS

BLOOD PRESSURE MANAGEMENT



5082-25

Adult Blood Pressure Cuff and Two-Tube Bladder



5082-200

Blood Pressure Connector (Tri-Purpose Type)



Two-Tube Blood Pressure Cuffs and Neoprene Bladders

For use with LARGE aneroid, wall aneroid, and mercury sphygmomanometers.

- Heavyweight nylon cuff
- Latex-free

5082-08 Infant Blood Pressure Cuff and Two-Tube Bladder

5082-24 Child Blood Pressure Cuff and Two-Tube Bladder

5082-25 Adult Blood Pressure Cuff and Two-Tube Bladder

5082-26 Large Adult Blood Pressure Cuff and Two-Tube Bladder

5082-78 Thigh Blood Pressure Cuff and Two-Tube Bladder

Two-Piece Inflation Systems

Completely assembled and ready to use.

- Heavyweight nylon material
- Latex-free
- Includes cuff, two-tube neoprene bladder, bulb and valve, and tri-purpose fitting

5082-22 Adult Inflation System with Two-Tube Bag

5082-23 Large Adult Inflation System with Two-Tube Bag

5082-45 Thigh Inflation System with Two-Tube Bag

5082-21 Child Inflation System with Two-Tube Bag

Two-Piece Inflation Systems: Accessories

5082-22 Adult Inflation System with Two-Tube Bag Two-Piece Inflation Systems Model Configuration

5082-23 Large Adult Inflation System with Two-Tube Bag Two-Piece Inflation Systems Model Configuration

5082-45 Thigh Inflation System with Two-Tube Bag Two-Piece Inflation Systems Model Configuration

5082-21 Child Inflation System with Two-Tube Bag

Two-Tube Neoprene Bladders

Latex-free neoprene for enhanced safety.

- Latex-free
- For use with pocket aneroid, wall aneroid, and mercury sphygmomanometers

5089-26 Black Neoprene Infant Inflation Bladder

5089-02 Black Neoprene Child Inflation Bladder

5089-01 Black Neoprene Adult Inflation Bladder

5089-22 Black Neoprene Large Adult Inflation Bladder

5089-19 Black Neoprene Thigh Inflation Bladder



5089-01
Black Neoprene Adult
Inflation Bladder

Inflation Bulbs and Valves

Durable blood pressure accessories in economical, individually wrapped packages.

Bulbs

5086-01 Small Inflation Bulb

5086-03 Large Inflation Bulb

5086-03H Large Inflation Bulbs (6 pack)

Bulb and Valve

5088-01 Small Inflation Bulb and Air-Release Valve

Air-Release Valve

5087-01 Blood Pressure Air-Release Valve



5088-01
Small Inflation Bulb and
Air-Release Valve



5086-03
Small Inflation Bulb



5086-01
Large Inflation Bulb



5087-01
Tycos® Air
Release Valve