



Electronic

Elite® Electronic Stethoscope	60
Elite Distributor	61
Elite Analyzer	61

Cardiology

Harvey™ Elite® Stethoscope	62
Harvey™ DLX™ Stethoscopes	64

Physical Diagnosis

Professional Adult Stethoscopes	66
Professional Pediatric Stethoscopes	67

General Exam Stethoscopes

Lightweight Double-Head Stethoscopes	68
--	----

Disposable Stethoscopes

Disposable Uniscopes	68
----------------------------	----



Elite® Electronic Stethoscope

With its excellent sound quality, the Elite helps you to hear more of the sounds you need to hear and less of what you don't. Using its two-position frequency selection switch, you can listen to either high or low frequencies to easily pick up low diastolic murmurs or high-pitched pulmonary sounds. In addition, accessories can be purchased that electronically capture, record and distribute body sounds or allow multiple users to co-listen, making the stethoscope perfect for teaching or consulting situations.

- Two-position filter switch lets you listen to high or low frequencies
- Adjustable binaurals and ear tips ensure all-day comfort
- Long-life lithium ion battery lasts up to 140 hours
- Latex-free
- Optional Distributor allows others to hear simultaneously
- Optional Analyzer stores and analyzes the sounds you hear right on your PC

5079-405 Elite Stethoscope

5079-405V Veterinary Elite Electronic Stethoscope (22" tube length)

Elite Electronic Stethoscope: Accessories

5079-426 Battery Cover for Elite Stethoscope

5079-427 Ear Tips for Elite Stethoscope

5079-425 3 V Lithium Battery CR123A for Elite Stethoscope

Elite® Distributor

The Elite Distributor allows you to connect up to five headsets to one stethoscope for group listening. This set-up is ideal for teaching, allowing all five individuals to hear a patient's body sounds at the same time. Each person can then make their own diagnosis based on the same sounds. It's a fast and simple way to train students and residents in auscultation and physical diagnosis.

5079-401 Elite Distributor (does not include headsets)

Elite Distributor: Accessories

5079-414 Headset (for use with the distributor)

5079-413 Distributor Cable



5079-414
Headset

Elite® Analyzer

The Elite Analyzer allows you to connect a stethoscope directly to a computer and record any exam. The software combines a phonocardiogram with a synchronized ECG signal, allowing the user to document any unusual findings or to monitor a patient's condition over time. Files can be emailed and/or stored in a patient's case file to be replaced at your convenience.

5079-430 Elite Analyzer (software with 3-Lead ECG)



5079-430
Elite Analyzer
(Computer not included)



5079-125
Harvey™ Elite® Stethoscope

Harvey™ Elite® Stethoscope

The Harvey Elite offers updated styling along with the superior comfort and auscultation capabilities that cardiologists have come to rely on. It is equipped with dual-bore, latex-free tubing that creates two uninterrupted sound channels and works with the optimally weighted, rugged, and stainless steel chestpieces. In addition, the Elite's diaphragm is designed to accentuate important sounds—all of this enables you to hear cardiac, pulmonary and vascular sounds.

- Dual lumen
- New longer 28" length
- Full-range cardiology/pulmonology stethoscope
- Double-head (bell and flat diaphragm) stainless steel chestpiece with nonchill rim
- Pediatric chestpieces available
- Rotatable stainless steel binaurals and interchangeable comfort sealing ear tips
- Triple-leaf binaural spring encased in polyurethane allows size adjustment while preventing breakage
- Choice of tubing length
- Latex-free

■	5079-122	Harvey Elite Stethoscope (22"), Black
■	5079-125	Harvey Elite Stethoscope (28"), Black
■	5079-270	Harvey Elite Stethoscope (28"), Burgundy
■	5079-271	Harvey Elite Stethoscope (28"), Navy
■	5079-284	Harvey Elite Stethoscope (28"), Forest Green
■	5079-125P	Harvey Elite Stethoscope (28"), Black, Pediatric Version

Harvey™ Elite® Stethoscope: Accessories

5079-113	Hard Ear Tips, Extra Large
5079-336	Comfort Sealing Ear Tips, Black
5079-337	Comfort Sealing Ear Tips (Large) Navy
5079-338	Comfort Sealing Ear Tips (Large), Burgundy
5079-339	Comfort Sealing Ear Tips(Large), Forest Green
5079-314	Comfort Sealing Ear Tips (Small), Gray
5079-233	Comfort Sealing Ear Tips (Large), Gray
5079-120	Tubing for 28" Model, Black
5079-272	Tubing for 28" Model, Burgundy
5079-305	Tubing for 28" Model, Forest Green
5079-273	Tubing for 28" Model, Navy
5079-274	Binaural/Spring Assembly, Burgundy
5079-121	Binaural/Spring Assembly, Black
5079-307	Binaural/Spring Assembly, Forest Green
5079-275	Binaural/Spring Assembly, Navy
5079-115	Pediatric Diaphragm Assembly, Black
5079-116	Adult Diaphragm Assembly, Black
5079-117	Pediatric Bell Nonchill Rim, Black
5079-118	Adult Bell Nonchill Rim, Black
5079-127	Pediatric Diaphragm Nonchill Rim, Black
5079-126	Adult Diaphragm Nonchill Rim, Black
5079-232	Accessory Kit, Black (Includes: One Pediatric Chestpiece, One Flat Diaphragm, One Corrugated Diaphragm, One Pair Soft Ear Tips, One Pair Extra Large Hard Ear Tips, One Pair Small Hard Ear Tips)
5079-282	Head Assembly, Black
5079-282G	Head Assembly, Gray
5079-104	Adult Flat Diaphragm Disc, Black
5079-124	Pediatric Diaphragm Disc Only
5079-05	Corrugated Diaphragm Disc Only
5079-119	Tubing for 22" Model, Black



5079-120
Black Tubing


Diaphragm Assemblies



5079-115
Pediatric



5079-116
Adult



5079-321
Harvey™ DLX,
Triple Head,
Black

Harvey™ DLX™ Stethoscopes

Both the double- and triple-head Harvey DLX stethoscopes are specifically designed to accentuate important sounds and enable you to hear cardiac, pulmonary, and vascular sounds without relying on difficult variable-pressure techniques. The Harvey DLX is equipped with dual-bore tubing that creates two uninterrupted sound channels and works with the optimally weighted brass chestpieces to transmit a full range of physical sounds. Providing clear, crisp sound, the Harvey DLX helps enhance your auscultation capabilities.

- New longer 28" length full-range, cardiology/pulmonary stethoscope
- Available with double-head (brass bell and flat diaphragm) or rotatable triple-head (brass bell and corrugated diaphragm) chestpieces
- Pediatric chestpieces available
- Chrome-plated brass chestpieces and binaurals
- Rotatable binaurals and interchangeable comfort sealing ear tips
- Triple-leaf binaural spring encased in polyurethane allows size adjustment while preventing breakage
- Latex-free
- Dual lumen



5079-325
Harvey™ DLX, Double
Head, Black

■	5079-321	Harvey DLX, Triple Head, Black (28")
■	5079-322	Harvey DLX, Triple Head, Burgundy (28")
■	5079-323	Harvey DLX, Triple Head, Navy (28")
■	5079-324	Harvey DLX, Triple Head, Forest Green (28")
■	5079-325	Harvey DLX, Double Head, Black (28")
■	5079-326	Harvey DLX, Double Head, Burgundy (28")
■	5079-327	Harvey DLX, Double Head, Navy (28")
■	5079-328	Harvey DLX, Double Head, Forest Green (28")
■	5079-325P	Harvey DLX, Double Head, Black, Pediatric (28")

FLAT DIAPHRAGM

Unique, tightly-contained diaphragm is optimized for higher-frequency, critically important heart sounds like murmurs, clicks, and ejection sounds.



CORRUGATED DIAPHRAGM

The corrugated diaphragm is especially good for detecting lower-frequency sounds, such as rumbles, gallops, and ejection sounds.



TRIPLE-HEAD MODELS

Get complete flexibility and efficiency by placing a bell, flat diaphragm, and corrugated diaphragm at your fingertips.



Harvey™ DLX™ Stethoscopes: Accessories

5079-336	Comfort Sealing Ear Tips (Large), Black
5079-337	Comfort Sealing Ear Tips (Large), Navy
5079-338	Comfort Sealing Ear Tips (Large), Burgundy
5079-339	Comfort Sealing Ear Tips (Large), Forest Green
5079-113	Hard Ear Tips, Extra Large
5079-314	Comfort Sealing Ear Tips (Small), Gray
5079-233	Comfort Sealing Ear Tips (Large), Gray
5079-108	Name Tag
5079-370	Diaphragm Assortment, Harvey Holographic, Adult, Four Colors
5079-272	Tubing for 28" Model, Burgundy
5079-120	Tubing for 28" Model, Black
5079-305	Tubing for 28" Model, Forest Green
5079-273	Tubing for 28" Model, Navy
5079-274	Binaural/Spring Assembly, Burgundy
5079-121	Binaural/Spring Assembly, Black
5079-307	Binaural/Spring Assembly, Forest Green
5079-275	Binaural/Spring Assembly, Navy
5079-104	Flat Diaphragm Only
5079-05	Corrugated Diaphragm Disc Only
5079-319	Diaphragm Assembly
5079-318	Bell Chestpiece Assembly
5079-316	Pediatric Chestpiece Assembly
5079-315	Harvey DLX Accessory Kit (Includes: Pediatric Chestpiece, Flat Diaphragm, Corrugated Diaphragm, One Pair of Soft Ear Tips, One Pair of Hard Ear Tips)
5079-124	Pediatric Diaphragm Disc Only



5079-370
Diaphragm Assortment



5079-120
Black Tubing








5079-135
Adult Stethoscope

Professional Adult Stethoscopes

With its minimal weight and custom-fitting headset and ear tips, the Professional is one of the most comfortable and versatile general-use stethoscopes manufactured today. Its acoustics allow for accurate auscultation of cardiac, pulmonary, and vascular sounds. And the chestpieces in the adult, pediatric, and anesthesiology versions of the Professional enable you to select a stethoscope that best meets your specific needs.

- General-use stethoscope
- Stainless steel double-head (bell and flat diaphragm) chestpiece with nonchill rim
- Rotatable binaurals and interchangeable comfort sealing ear tips
- Latex-free

	5079-135 Adult Stethoscope, 28" (71 cm), Black
	5079-137 Adult Stethoscope, 28" (71 cm), Navy
	5079-289 Adult Stethoscope, 28" (71 cm), Blue
	5079-139 Adult Stethoscope, 28" (71 cm), Burgundy
	5079-285 Adult Stethoscope, 28" (71 cm), Forest Green

Professional Adult Stethoscope: Accessories





5079-123	Adult Diaphragm Disc Only
5079-124	Pediatric Diaphragm Disc Only
5079-180	Adult Bell Nonchill Ring, Black
5079-181	Adult Bell Nonchill Ring, Gray
5079-182	Pediatric Bell Nonchill Ring, Black
5079-183	Pediatric Bell Nonchill Ring, Gray
5079-126	Adult Diaphragm Nonchill Ring, Black
5079-184	Adult Diaphragm Nonchill Ring, Gray
5079-127	Pediatric Diaphragm Nonchill Ring, Black
5079-185	Pediatric Diaphragm Nonchill Ring, Gray
5079-195	Binaural/Spring Assembly and Tubing, 28", Black
5079-196	Binaural/Spring Assembly and Tubing, 28", Green
5079-197	Binaural/Spring Assembly and Tubing, 28", Navy
5079-199	Binaural/Spring Assembly and Tubing, 28", Burgundy

5079-210	Double Head Chestpiece Assembly (Adult), Gray
5079-211	Double Head Chestpiece Assembly (Adult), Black
5079-214	Double Head Chestpiece Assembly (Pediatric), Black
5079-215	Double Head Chestpiece Assembly (Pediatric), Gray
5079-320	Double Head Chestpiece Assembly (Infant), Black
5079-172	Hard Ear Tips (Medium), Pair
5079-173	Hard Ear Tips (Large), Pair
5079-313	Comfort Sealing Ear Tips (Small), Gray
5079-366	Comfort Sealing Ear Tips (Large), Black
5079-170	Comfort Sealing Ear Tips (Large), Gray
5079-367	Comfort Sealing Ear Tips (Large), Navy
5079-368	Comfort Sealing Ear Tips (Large), Burgundy
5079-369	Comfort Sealing Ear Tips (Large), Forest Green

Professional Pediatric Stethoscopes

With its minimal weight and custom-fitting headset and ear tips, the Professional is one of the most comfortable and versatile general-use stethoscopes manufactured today. Its acoustics allow for accurate auscultation of cardiac, pulmonary, and vascular sounds.

- Meets all diagnostic needs, from detecting Korotkoff sounds of blood pressure to analyzing subtle heart murmurs
- Stainless steel chestpiece intensifies body sounds transmitted through tissues
- Comfort sealing ear tips depress and conform to the shape of ear canals for soft yet airtight fit
- Latex-free

5079-145	Pediatric Stethoscope, 28" (71 cm), Black	
5079-147	Pediatric Stethoscope, 28" (71 cm), Navy	
5079-149	Pediatric Stethoscope, 28" (71 cm), Burgundy	
5079-287	Pediatric Stethoscope, 28" (71 cm), Forest Green	



5079-145
Pediatric Stethoscope



5079-75,74,73
Lightweight Double-Head
Stethoscopes

Lightweight Double-Head Stethoscopes

When weight and economy are of primary concern, the Lightweight is the stethoscope to choose. At 4 oz, the Lightweight is one of the lightest stethoscopes available. It rests comfortably in your lab coat or draped around your neck, and offers the superior auscultation capabilities that are the hallmark of Welch Allyn Stethoscopes.

- General-use stethoscope
- Lightweight chestpieces with nonchill rims
- Adjustable chrome-plated binaurals
- Latex-free, flexible, single-channel tubing

	5079-73	Lightweight Double-Head Stethoscope, Velvet Black
	5079-76	Lightweight Double-Head Stethoscope, Misty Blue
	5079-74	Lightweight Double-Head Stethoscope, Poppy Red
	5079-75	Lightweight Double-Head Stethoscope, Dawn Gray



17461
Disposable Uniscope, yellow

Disposable Uniscopes

Designed as a single-use stethoscope for situations where contamination issues are of concern, the Welch Allyn Uniscope disposable stethoscope enables you to make accurate diagnoses without compromising your primary stethoscope.

- General-use stethoscope
- For use in situations where universal precautions are mandated
- Single-head chestpiece
- Pediatric and anesthesia models available
- Latex-free tubing
- MRI compatible

	17461	Disposable Uniscope, Yellow
	17462P	Disposable Uniscope, Pediatric, Light Blue

Welch Allyn—Stethoscope Specifications Chart

	ELITE ELECTRONIC STETHOSCOPE	CARDIOLOGY GRADE		PROFESSIONAL GRADE			DISPOSABLE
PRODUCT DESCRIPTION	ELITE ELECTRONIC	HARVEY DLX (28")	HARVEY ELITE DOUBLE HEAD (28") (22")	PROFESSIONAL ADULT	PROFESSIONAL PEDIATRIC	LIGHTWEIGHT	UNISCOPE
PRODUCT NUMBERS	5079-405 5079-430 (Analyzer) 5079-401 (Distributor)	5079-321 Triple, Black 5079-322 Triple, Burgundy 5079-323 Triple, Navy 5079-324 Triple, Green 5079-325 Double, Black 5079-326 Double, Burgundy 5079-327 Double, Navy 5079-328 Double, Green	5079-122 Elite, Black (22") 5079-125 Elite, Black 5079-270 Elite, Burgundy 5079-271 Elite, Navy 5079-284 Elite, Green 5079-125P Elite, Ped Blk	5079-135 Black 5079-137 Navy 5079-289 Blue 5079-139 Burgundy 5079-285 Forest	5079-145 Black 5079-147 Navy 5079-149 Burgundy 5079-287 Forest	5079-73 Velvet Black 5079-74 Poppy Red 5079-75 Dawn Gray 5079-76 Misty Blue	17461 Adult, Yellow 17462P Pediatric, Blue
PRODUCT TYPE	Electronic Stethoscope	Double/Triple Head	Double Head	Double Head	Double Head	Double Head	Single Head
HEAD DIAMETER INCHES (mm)	1.5 (38.5)	1.9 (48.3)	1.9 (48.3)	1.75 (45)	1.2 (30)	1.9 (48.3)	1.9 (49)
OVERALL LENGTH INCHES (mm)	32.5 (825.5)	28 (711)	28 (711)	28 (711)	28 (711)	32 (812.5)	29.5 (750) 27 (686)
TUBING LENGTH INCHES (mm)	21.5 (546.1)	17 (431.8)	17 (431.8)	21.5 (546.1)	21.5 (546.1)	24 (609.6)	19 (483)
LATEX FREE (Y/N)	Yes	Yes	Yes	Yes	Yes	Yes	Yes
CHESTPIECE MATERIAL	Plastic	Brass	Stainless Steel	Stainless Steel	Stainless Steel	Chrome Plated Zinc	ABS
DIAPHRAGM MATERIAL	Plastic	Polycarbonate	Polycarbonate	Clear Polycarbonate	Iten Resiten G-10/FR4T Glass Epoxy Sheet	PVC	Polypropylene
DIAPHRAGM STYLE	Flat	Flat or Corrugated	Flat or Corrugated	Flat	Flat	Flat	Flat
Adjustable Binaural (Y/N)	Yes	Yes	Yes	Yes	Yes	Yes	No
EAR TIPS (SOFT/HARD)	Soft	Soft/Hard	Soft/Hard	Soft/Hard	Soft/Hard	Hard	Hard
TYPE OF TUBING (SINGLE/DUAL LUMEN)	Single Lumen	Dual Lumen	Dual Lumen	Single Lumen	Single Lumen	Single Lumen	Single Lumen
NONCHILL RING (Y/N)	No	Yes (Brass)	Yes	Yes	Yes	Yes	No
SHIPS ASSEMBLED (Y/N)	Yes	Yes	Yes	Yes	Yes	No	Yes
PACKAGING	Printed Box	Printed Box	Printed Box	Printed Box	Printed Box	Plain White Box 1-piece	Plastic Bag
DESCRIPTION	Meets all diagnostic needs of the demanding professional. Provides the user with the ability to filter and amplify heart and lung sounds during auscultation. External 2.5 mm jack can be used to connect the stethoscope to an external device such as the Analyzer or the Distributor. Excellent tool for those experiencing hearing loss. Uses one 3 V lithium ion battery. Latex-free.	Designed to meet the demanding needs of cardiologists, this stethoscope has fully adjustable binaurals and triple leaf spring for comfort. The head is made of trumpet brass because of its ability to transmit sounds accurately. The dual-lumen tubing is latex-free.		Meets all diagnostic needs, from detecting Korotkoff sounds of blood pressure to analyzing subtle heart murmurs. Stainless-steel chestpiece intensifies body sounds transmitted through tissues. Comfort sealing ear tips depress and conform to the shape of ear canals for soft yet airtight fit. Latex-free.		Quality design for blood pressure measurement and general auscultation. Lightweight and flexible single-channel tubing provides comfortable use over extended periods. Adjustable chrome-plated binaurals. Latex-free.	The Uniscope is a single-patient-use stethoscope for situations where cross-contamination is a concern.