





### Electronic Thermometry

SureTemp® Plus 690.....	134
SureTemp Plus 692 .....	135

### Ear Thermometry

Braun® ThermoScan® PRO 4000.....	137
----------------------------------	-----



**01690-200**  
SureTemp® Plus 690  
Thermometer



**01690-400**  
SureTemp® Plus 690 Electronic  
Thermometer, Wall Mount, Oral  
Probe with Oral Probe Well

### SureTemp® Plus 690 Electronic Thermometer

Easy to use with intuitive operation, reduces the risk of cross-contamination with removable probe well and provides more reliable readings with improved accuracy. Takes approximately 4-to-6 second oral temperatures, 10-to-13 second pediatric axillary and rectal temperatures and 12-to-15 second adult axillary temperatures.

- Ergonomic design—comfortable to hold, easy to use, and very durable
- Interchangeable, removable probe wells reduce the risk of cross-contamination
- Large LCD displays temperatures in Fahrenheit or Celsius
- Last temperature recall
- Convenient storage for 25 probe covers
- Takes approximately 6,000 readings on 3 AA batteries

**01690-200** SureTemp Plus 690 Electronic Thermometer, 4 ft Cord and Oral Probe with Oral Probe Well

**01690-201** SureTemp Plus 690 Electronic Thermometer, 4 ft Cord and Rectal Probe with Rectal Probe Well

**01690-300** SureTemp Plus 690 Electronic Thermometer, Wall Mount, 9 ft Cord and Oral Probe with Oral Probe Well

**01690-301** SureTemp Plus 690 Electronic Thermometer, Wall Mount, 9 ft Cord and Rectal Probe with Rectal Probe Well

**01690-400** SureTemp Plus 690 Electronic Thermometer, Wall Mount, 4 ft Cord and Oral Probe with Oral Probe Well

**01690-700** SureTemp Plus 690 Electronic Thermometer, Wall Mount, 9 ft Cord and Oral Probe with Oral Probe Well, with Rolling Stand and Holder for Spare Probe

## SureTemp® Plus 692 Electronic Thermometer

Easy to use with intuitive operation, reduces the risk of cross-contamination with removable probe well and provides more reliable readings with improved accuracy. Takes approximately 4-to-6 second oral temperatures, 10-to-13 second pediatric axillary and rectal temperatures and 12-to-15 second adult axillary temperatures.

- Ergonomic design—comfortable to hold, easy to use, and very durable
- Interchangeable, removable probe wells reduce the risk of cross-contamination
- Large back lit LCD display with location identification capability
- Includes a 60-second pulse timer for counting pulses and respirations
- Mechanical or electronic security options to help prevent theft
- Takes approximately 5,000 readings on 3 “AA” batteries
- Wall mount standard with the SureTemp Plus 692 model



SureTemp® Plus 692  
Electronic Thermometer

- |                  |  |
|------------------|--|
| <b>01692-200</b> | SureTemp Plus 692 Electronic Thermometer with Wall Mount, Security System with ID Location Field, 4 ft Cord and Oral Probe with Oral Probe Well  |
| <b>01692-201</b> | SureTemp Plus 692 Electronic Thermometer with Wall Mount, Security System with ID Location Field, 4 ft Cord and Rectal Probe with Rectal Probe Well  |
| <b>01692-300</b> | SureTemp Plus 692 Electronic Thermometer, with Wall Mount, Security System with ID location field, 9 ft Cord and Oral Probe with Oral Probe Well   |
| <b>01692-301</b> | SureTemp Plus 692 Electronic Thermometer, with Wall Mount, Security System with ID location field, 9 ft Cord and Rectal Probe with Rectal Probe Well   |
| <b>01692-700</b> | SureTemp Plus 692 Electronic Thermometer with Wall Mount, Security System with ID location field, 9 ft Oral Probe with Oral Probe Well, with Rolling Stand and Holder for Spare Probe and Well |



**05031-101**  
Disposable Probe Covers



**21330-0000**  
Holder for Spare Probe and Well



**02893-000**  
Oral Probe and Well Kit, 4 ft



**02891-1000**  
Probe Well, Red,  
Rectal

### SureTemp Plus: Accessories

**05031-101** Disposable Probe Covers (1,000 covers, packaged 25 per box)

**05031-110** Disposable Probe Covers (10,000 covers, packaged 25 per box)

**21326-0000** Wall Holder

**21330-0000** Holder for Spare Probe and Probe Well

**06138-000** Calibration Key

**02892-000** Probe and Well Kit, 4 ft Rectal

**02892-100** Probe and Well Kit, 9 ft Rectal

**02893-000** Probe and Well Kit, 4 ft Oral/Axillary

**02893-100** Probe and Well Kit, 9 ft Oral/Axillary

**02891-0000** Probe Well, Blue, for Oral/Axillary

**02891-1000** Probe Well, Red, for Rectal

**08273-000** Thermometer Roll Stand with Holder for Spare Probe and Well

**21333-0000** Bracket, M600 Stand Upgrade with Holder for Spare Probe and Well

**01802-110** 9600 Plus Calibration Tester

### Braun® ThermoScan® PRO 4000 Ear Thermometer

Provides the speed and accuracy you need for all types of patients. The exclusive ExacTemp™ technology provides an active user feedback system to ensure proper positioning and improved accuracy.

- Improved accuracy and repeatability through a preheated probe tip and sensor
- Probe cover detection system ensures precise infrared heat transfer for accurate temperatures
- Automatic probe cover eject button for quick and easy disposal of used probe covers
- Easy-to-read LCD displays temperatures in Fahrenheit or Celsius

**04000-200** Braun ThermoScan PRO 4000 Ear Thermometer (with 2 "AA" batteries)

**04000-300** Braun ThermoScan PRO 4000 Ear Thermometer with Recharging Base Station

**04000-400** Braun ThermoScan PRO 4000 Ear Thermometer with Recharging Base Station and Rolling Thermometer Stand

### Braun® ThermoScan PRO 4000: Accessories

**05075-800** Braun ThermoScan PRO 4000 Disposable Probe Covers (800 covers, packaged 20 per box)

**05075-005** Braun ThermoScan PRO 4000 Disposable Probe Covers (5,000 covers, packaged 20 per box)

**24001-0000** Rolling Thermometer Stand and Recharging Base Station

**24001-1000** Recharging Base Station and Probe Cover Refill Dispenser

**53020-0000** Rechargeable Battery Pack

**01802-110** 9600 Plus Calibration Tester

**701640** Replacement Battery Cover

**713420** Replacement Cradle

**80199-0000** Replacement Power Cord for PRO 4000 Base Station



**04000-200**  
Braun® ThermoScan® PRO  
4000 Ear Thermometer

**04000-400**  
Braun® ThermoScan PRO 4000  
Ear Thermometer with Rolling  
Thermometer Stand



**04000-300**  
Braun® ThermoScan® PRO  
4000 Ear Thermometer with  
Recharging Base Station



**05075-800**  
Disposable Probe Cover



**53020-0000**  
Rechargeable Battery Pack



**24001-1000**  
Recharging Base Station