

Spot Vital Signs® Monitor Competency Checklist

Name (print): _____ Unit: _____ Date: _____

- _____ Locate the Power button and turn the Spot device on.
- _____ Obtain a NIBP (non-invasive blood pressure) measurement.
- _____ True or False. If a patient falls into the range between a regular and large Adult blood pressure cuff, you should use the larger.
- _____ True or False. The Spot Vital Signs® with SureTemp thermometry can be used to take oral, axillary, or rectal temperatures.
- _____ Switch the thermometer mode to oral. What is the proper probe positioning in the mouth?
- _____ True or False. Placing the pulse oximeter sensor on the same limb as the blood pressure cuff will have no effect on the SpO2 reading.
- _____ What do the letters 'Z Z Z' across the top of the monitor display mean? Bring the Spot Vital Signs® out of standby mode.
- _____ If the device screen is blank, how do you recall the last vital signs measurement?
- _____ True or False. You should press Clear (C) before taking vital signs on the next patient.